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Algeness® is an absorbable filler indicated for the correction of tissue and / or hypodermic atrophies. **Algeness®** is an injectable gel consisting of agarose.

Algeness® is contained in a special "total emptying" double syringe that allows you to mix the product before use and to use all the gel contained in the syringe.

COMPOSITION

Algeness® LD contains agarose (1%), water (89,1%) and phosphate buffer (9,9%).

Algeness® HD contains agarose (1,5%), water (88,7%) and phosphate buffer (9,8%).

Algeness® VL contains agarose (2,5%), non cross-linked hyaluronic acid (0,5%), water (87,3%) and phosphate buffer (9,7%).

Algeness® DF contains agarose (3,5%), non cross-linked hyaluronic acid (0,4%), water (86,4%) and phosphate buffer (9,7%).

Agarose, chosen for its physical and rheological properties suitable for the application site, is able to form an extrudable gel consisting of a three-dimensional matrix capable of mechanically retaining molecules and organic liquids in dynamic equilibrium with its absorption. Therefore, once the agarose has been injected, thanks to its structure, it is able to give viscoelastic support to the extracellular matrix of the tissue by mechanically filling the spaces involved in subdermal atrophy.

Algeness®, injected into the subdermal tissue, provides a viscoelastic supplement to the matrix. Agarose restores the lost tissue volume of both the adipose tissue and the connective stroma, ensuring biocompatibility with the extracellular matrix and the harmony of natural forms.

INTENDED USE

Algeness® is effective for the treatment of cases of dermal atrophy and connective tissue deficits. **Algeness®** can also be used for deep filling of cutaneous tissue.

Treatment with **Algeness** is part of reconstructive surgery which aims to correct congenital malformations or functions compromised by trauma, disease, accidents, burns, wounds, often accompanied by loss of substance.

For **Algeness**, the clinical advantage is the reduction of the psychological impact, which can have pathological consequences, caused by the morphological alterations of the patient's body.

Algeness® LD (low density - subdermal) is indicated for repairs and superficial subdermal corrections surface repairs and subdermal corrections of hypotrophy or tissue lesions such as in the perioral area.

Algeness® HD (medium density - subdermal) is indicated for repairs and superficial subdermal corrections due to hypotrophy or tissue lesions in the perioral area.

Algeness® VL (medium / high density - deep subdermal) is indicated in the volumetric restoration for all forms of deep tissue hypotrophy with severe or modest loss of tone of the deep ligaments such as in the suborbital areas and in the anterior mandibular area.

Algeness® DF (high density - deep subdermal) is indicated for all forms of volumetric restorations for all forms of severe deep tissue hypotrophy with severe loss of volume and tone of the deep ligaments such as in the suborbital areas and in the anterior mandibular area.

INDICATIONS AND USE

Prior to use, mix the gel back and forth between the two syringes at least ten times (10) for improving consistency and better results. Then, disconnect the empty syringe with the connector and attach the supplied needle to the pre-filled syringe.

Algeness® is a slowly absorbed filler and is totally biocompatible, it is indicated for the treatment of atrophic and hyperkeratotic changes of the skin and for the correction of wrinkles, skin folds, scars, or for problems associated with traumatic injuries.

The treatment must only be performed by competent doctors trained in the use of injection techniques.

Before using the product, the patient must be informed about the indications, contraindications and possible side effects that may occur with use.

Before the injection, it is necessary to carry out a thorough cleansing and complete disinfection of the implant site. It is recommended to maintain a correct asepsis of the treated area before, during and after the insertion of the implant as well as to ensure a proper environment (clinic).

Position **Algeness®** in the medium-deep subdermis under the site to be corrected or deeper, depending on the concentration of the agarose. The quantity to be injected is left to the judgment of the doctor.

Algeness® should not be used in excess (overcorrection).

The correction can be performed with the usual techniques: linear threading technique, deep linear or fan technique. The injection must be carried out slowly for better positioning of the product in the desired areas with less trauma to the tissues.

The product is extruded by applying continuous and constant pressure on the plunger of the syringe while withdrawing the needle.

For an optimal correction, after implantation, always massage the treated area to make a homogeneous distribution of the injected material.

It is recommended to apply a cold compress to the treatment area to reduce any local undesirable reactions.

For proper maintenance of the result, remind the patient of the importance of a clinical follow-up a few months after treatment.

The patient must also be informed that the treatment can be suspended at any time at their request.

Inform the patient they must not apply cosmetics to the implant site for the next twelve hours and not to expose themselves to direct heat for the next few days (for example: exposure to the sun or to UVA and UVB rays, use of hair dryers or of hair dryers, heat radiated by fireplaces, saunas, etc.)

CONTRAINDICATIONS

The use of **Algeness®** is contraindicated in all cases other than those specified in the INDICATIONS of the product, therefore, it must not be injected superficially at the level of the dermis or as a bolus.

Prior to the injection, assure that the patient does not experience particular hypersensitivity to any of the components in the product.

As with all fillers, the use of the product for the correction of particularly vascularized areas is not recommended, as it can increase the risk of compression and occlusion of the vessels and related phenomenology.

Algeness® must not be injected into blood vessels, as it could cause the onset of occlusion, local tissue necrosis or emboli.

Algeness® must not be injected in subjects with acute or chronic skin diseases in the areas to be corrected or in the immediate vicinity and, for precautionary reasons, in subjects with a history of anaphylactic reactions or severe allergies, in patients with severe organ or systemic diseases, including autoimmune disorders. The product is not intended for use on children, pregnant and or lactating women. It is recommended to avoid the combined use with other substances such as crosslinked fillers in the same treatment area.

PRECAUTIONS

As with all percutaneous procedures, the implant can entail risks of infection if the product is injected into anatomical areas where inflammatory or infectious processes are present, or without proper cleansing and disinfection of the area to be corrected.

In patients with bleeding and / or coagulation disorders or during treatments with anticoagulants, the product should be used with caution, as the injection may frequently cause local bleeding or bruising.

ADVERSE REACTIONS

Following the implantation, edematous reactions or skin redness, of a mild or moderate degree, rarely occur, however these will resolve completely within a few days.

Even if the application of **Algeness®** is not painful, the possibility of rare reports of transient pain due to the injection cannot be excluded, especially in particularly sensitive areas.

WARNINGS

If you do not perform a thorough massage following the injection, it is possible the appearance of lumps in the treated area.

PRODUCT IS FOR SINGLE USE ONLY. STERILIZED product. Product must not be re-sterilized.

If the stickers on the inner package, which indicate the sterility of the product are NOT red, DO NOT use the product. In case of any irregularities with the packaging, isolate the non-compliant packaging and dispose of as special waste. Never use the product if it has not been stored correctly or if it has been frozen.

The unused material left in the syringe must be eliminated (disposed of) after treatment to avoid any risk of cross-infection due to use on other patients.

The doctor and / or the patient must report any serious incidents occurring in relation to the device to the manufacturer Ghimas and to the competent authority of the Member States where the user and / or patient is established.

FOR PROFESSIONAL USE ONLY.

DO NOT USE the product after the expiry date shown on the package.

PACKAGING

The product is supplied in a pouch with a special double syringe Luer Lock "total emptying". Each box contains two pouches and four sterile needles.

The volume contained in the syringe is indicated on the pouch label.

STORAGE

Conserve the product at room temperature, between 6°C and 30°C. If necessary, the product may be stored in a refrigerator, the product must not be frozen.

On the basis of MDR 2017/745 and subsequent modifications, the physician must ensure the traceability of the product, i.e. he must be able to trace the patient to whom he applied the product. It is recommended to attach the product label to the patient's medical record. Consult the attachment of the instructions for completing the implant card, take the inserted card and deliver it completed to the patient.

It is recommended to carefully read the Summary relating to safety and clinical performance of **Algeness** published on the Ghimas website, www.ghimas.it, in the dedicated area and/or on the European database "Eudamed".

The product must be used according to the instructions for use.

The manufacturer assumes no responsibility for any damages deriving from different and inappropriate use. The user is required to personally check the suitability of the product for the application he intends to use, especially if they are not listed in the instructions for use.

These instructions were drafted on: 23/06/2023

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Advanced Aesthetic Technologies Inc. – Brookline MA, USA